


## Live click labelling with TriPPPPro

For any questions regarding this protocol, feel free to send an email to [info@deep-dv.org](mailto:info@deep-dv.org)

### Purpose

The TriPPPPro system facilitates live imaging of cellular and viral genomes through the use of a pronucleotide linked to a dienophile and a fluorogenic linked to a diene. These components undergo a catalysator-free inverse electron demand Diels Alder reaction with rapid reaction kinetics without the need for fixation and permeabilization of the cells, allowing the performance of live cell experiments (Sterrenberg, Stalling et al., *Angewandte Chemie* 2023, <https://doi.org/10.1002/ange.202308271>).



**Figure 1:** Graphical Scheme of labeling genomes with the TriPPPPro system.

### Notes

- To obtain the TriPPPPro and their corresponding dyes, please reach out to Iven Knaack (Research Group of Prof. Dr. Chris Meier, Institute of Organic Chemistry, University of Hamburg).
- Cells should have a confluency of 70-80% before proceeding with the click labeling.
- Samples can be fixed in 4% PFA after completing the click reaction.

### Reagents

Reagent/Material	Description	Reference/Vendor	Cat.no.
<b>μ-slide 8 well high glass bottom</b>	Glass coverslip bottom microscopy 8-well chamber slide	ibidi	80807
<b>DMEM</b>	Supplement with 10% FBS	Sigma-Aldrich	D6429
<b>FBS Superior</b>		Merck	F7524
<b>VP SFM (1x)</b>	Virus production serum-free medium, supplement with 20 ml GlutaMAX™ Supplement	ThermoFisher	11681020
<b>GlutaMAX™ Supplement</b>		ThermoFisher	35050061
<b>Hoechst 33342</b>	Stains DNA	ThermoFisher	H3569
<b>TriPPPPro 2TCOa-dCTP</b>	1*10 <sup>-2</sup> mol/l	UHH, AK Meier	
<b>Click-Fluorophore</b>	1.73*10 <sup>-3</sup> mol/l, HD656, excitation maximum: 656 nm	UHH, AK Meier	IK341

# Live click labelling with TriPPPPro Protocol

## Day 0

Seed cells (e.g. Vero cells,  $2,5 \cdot 10^4$  cells / well) in 10% FBS DMEM in ibidi 8-well plate

## Day 1


-x h Infect cells if applicable

0 h Wash cells 3x in serum-free medium (VP SFM)  
Add 5  $\mu$ M TriPPPPro in serum-free medium (VP SFM)  
-> 1000  $\mu$ l VP SFM + 0.5  $\mu$ l TriPPPPro

2 h Wash cells 3x in in serum-free medium (VP SFM)  
Add 5  $\mu$ M IK341 dye and Hoechst (1:1000, if necessary) in serum-free medium (VP SFM)  
-> 1000  $\mu$ l VP SFM + 2.9  $\mu$ l IK341

2.30 h Live cell imaging with a confocal microscope, e.g. Spinning Disk microscope  
Use a 647 nm laser for imaging the TriPPPPro-labelled genomes

## Results



**Figure 2: Live click labelling of nascent viral genomes in living cells.** Vero cells were infected with HSV-1 VP26\_mCherry at an MOI of 10. 1 hpi medium was changed and 4 hpi cells were pulsed with 5  $\mu$ M TriPPPPro-dCTP for 2 hours. Scalebar 10  $\mu$ m.